


Instrument Overview

- Hardware setup can be divided into a few key assemblies:
 - Microscope Assembly
 - CE Detection Assembly
 - CE Outlet Assembly




CE Detection Assembly

Custom Parts



Custom Material List

| Object | Material | Fabrication |
|------------------------------|----------|--------------|
| SiPM Mount | PLA/PETG | 3D Print |
| Capillary Arm | PLA/PETG | 3D Print |
| Cantilever Support | Aluminum | Mill |
| Optical Mount | PLA/PETG | 3D Print |
| Smooth Post Mount | PLA/PETG | 3D Print |
| Capillary Alignment Assembly | PLA/PETG | 3D Print |
| RGB-PCB | FR4 | EDA/PCB Fab. |



CE Detection Assembly Purchased Components


| | Component | Company | Part number | Notes |
|---|--|---------------|-------------------------|--------------------|
| A | 40X Objective 0.75 NA | Nikon | | |
| B | 4-Way Mounting Cage | Thorlabs | C4W | |
| C | Cage Dichroic Filter Mount | Thorlabs | FFM1 | |
| D | Rotatable Cover Plate | Thorlabs | B3C | |
| E | Blank Cover Plate | Thorlabs | B1C | |
| F | 10 mm Lens Tube | Thorlabs | SM1L10 | |
| G | External SM1 threads to internal M25x0.75 Thread (objective) | Thorlabs | SM1A12 | |
| H | External C-mount to SM1 thread | Thorlabs | SM1A39 | *PMT |
| I | Calibrated Iris Diaphragm | Thorlabs | SM1D12C | |
| J | SM1 Threaded Manual Shutter | Thorlabs | SM1SH1 | *PMT |
| K | End Cap with SM1 Threads | Thorlabs | SM1CP2 | |
| L | Non Roating adjustable SM1 Lens Tube | Thorlabs | SM1ZM | |
| M | Dichroic (488 nm) | Chroma | ZT488rdc-UF2 (IN036384) | 25.5 x 36 x 2mm |
| N | Bandpass Filter (525 +/- 25 nm) | Chroma | ET525/50m (IN005750) | 25 mm dia mounted |
| O | Notch Filter (488 nm) | Chroma | ZET488NF | 25 mm dia mounted |
| Q | Fiber Receptacle Collimators | OZ Optics | HPUCO | |
| P | gasket tape | McMaster-Carr | 95705K133 | *PMT .02" thick |
| R | Optical post | Generic | | 25 mm Diameter |



* Items needed for Hamamatsu PMT detector. Can be ignored when using SiPM

CE Outlet Assembly Custom Parts

| Object | Material | Fabrication |
|---------------------------|------------|----------------|
| Outlet Vial Arm | PLA/PETG | 3D Print |
| Outlet Vial Clamp | PLA/PETG | 3D Print |
| Lead Screw Smooth Bracket | PLA/PETG | 3D Print |
| Motor Smooth Bracket | PLA/PETG | 3D Print |
| Optical Post Bracket | PLA/PETG | 3D Print |
| Limit Switch | Breadboard | Snips/Scissors |




CE Outlet Assembly

Purchased Components

| Object | P/N | Supplier | Notes |
|-----------------------------|------------------|------------------|---|
| 2.54 mm Pitch Connectors | HDB-635 | Hilitchi, Amazon | Any standard crimped 2.54 mm connectors |
| Limit Switch | DBWDKG-FT01 | MUZHI, Amazon | Any standard SPDT Switch |
| Lead Screw w/ Motor and nut | RDX3D060 | Redrex, Amazon | Any leadscrew + motor combo* |
| Vial, Cap, and Septum | 27340 | MilliporeSigma | Can be swapped with other 4 mL open top screw cap vials with PTFE/Silicon septum. |
| M3 Nylon Screw | 92492A721 | McMasterCarr | 16 mm L |
| M3 Screws + Hex Nuts | M3-0420F-SS305 | Iexcel, Amazon | Generic assortment of lengths |
| Arduino Uno | | Arduino | |
| PowerSTEP01 | X-NUCLEO-IHM03A1 | Mouser | |
| 12 V Power Supply | | SmoTecQ, Amazon | The driver can power 12-48V so choose a power supply suitable for your motor. |


* 8 mm lead on the lead screw and 1.8° step for the motor will not require adjusting the OutletControl.ino code.




Exploded view of the outlet assembly.

Limit Switch for Outlet

- Simple through hole bread board can be used.
- Most generic single pull double throw (SPDT) limit switches will work.
- Breadboard can be cut with box cutters or snips. To rough dimensions specified.
- Breadboard will need to have two holes drilled for 3M mounting spaced 35 mm apart, 3.1mm diameter.




Wiring diagram for Limit Switch



Breadboard dimensions in mm

Outlet Arduino Pin Assignments

- Place the X-Nucleo ontop of the Arduino Uno.
 - Pin spacing is the same for Nucleo and Arduino Uno.*
- 3 Wires are reserved for the limit switch (A4, A5, D3).
- 3 Wires are reserved for the inlet RGB LED (A1, A2, A3)
- 3 Wires are reserved for the Solenoid (D4, D6, & D7)
- Wire the stepper motor to the Nucleo output screw terminals (A+/A-, B+/B-).
 - Each coil of the stepper motor is indicated by A or B. Identify which leads correspond to the same coil using an Ohm meter.
- Connect the 12 V power supply to the Nucleo power screw terminals (Labeled GND and +V).



<https://store.arduino.cc/usa/arduino-uno-rev3>